

leXsolar-BioFuel Large



Student's Manual

Layout diagram leXsolar-BioFuel Large

Item-No.1702

Bestückungsplan leXsolar-BioFuel Large

Art.-Nr.1702



- 1 L2-06-079 Areometer
L2-06-079 Dichtemessspindel
- 2 3xL2-06-085 Pasteur pipette
3xL2-06-085 Pasteurpipette
- 3 L2-06-072 Alcoholmeter
L2-06-072 Alkoholmeter
- 4 L2-06-070 Distilling head
L2-06-070 Destillieraufsatz

- 5 L2-06-076 Airlock
L2-06-076 Gärspund
- 6 L2-06-084 Grip stopper
L2-06-084 Lamellen Stopfen
- 7 L2-06-087 Syringe 2 ml
L2-06-087 Spritze 2 ml
- 8 3xL2-06-083 Test tubes
3xL2-06-083 Reagenzgläser

Version number
Versionsnummer


L3-03-142_06.05.2015

Layout diagram leXsolar-BioFuel Large

Item-No.1702

Bestückungsplan leXsolar-BioFuel Large

Art.-Nr.1702



- | | | | |
|----|--|----|--|
| 9 | L2-06-016 Laboratory thermometer
L2-06-016 Laborthermometer | 14 | L2-06-075 Erlenmeyer flask 1000 ml
L2-06-075 Erlenmeyerkolben 1000 ml |
| 10 | L2-06-071 Condenser
L2-06-071 Kühler | 15 | L2-06-082 Beaker 250 ml
L2-06-082 Becherglas 250 ml |
| 11 | L2-06-086 Measuring cylinder 100 ml
L2-06-086 Messzylinder 100 ml | 16 | 1702-02 Yeast
1702-02 Hefe |
| 12 | 1100-27 Motor module without gear
with L2-02-017 Propeller
1100-27 Motormodul ohne Getriebe
mit L2-02-017 Propeller | 17 | 1702-01 Plug with hose
1702-01 Stopfen mit Schlauch |
| 13 | 1700-01 Ethanol fuel cell module
1700-01 Ethanol-Brennstoffzellenmodul | 18 | 1100-23 Potentiometer module
1100-23 Potentiometermodul |
| | | 19 | L3-03-016 leXsolar-CD
L3-03-016 leXsolar-CD |

leXsolar-Bio Fuel Large

Student`s Manual

Content

1.1 Production of Biodiesel (FAME) from fats & oils	4
1.2 Extraction of fats from groceries or oil plants	7
1.3 Differences between cooking oil and FAME	9
2.1 Fermentation on the basis of different sugars.....	10
2.2 Proof of CO ₂ by a calcium hydroxide solution	12
2.3 Link between reaction speed and temperature	14
3 Distillation of fermented mash	16
4.1 Working principle of the fuel cell	22
4.2 Recording the characteristic curve of a 15 % solution	23
4.3 Temperature dependence of the power	25
4.4 Concentration dependence of the power and the characteristic curve	27




1.1 Production of Biodiesel (FAME) from fats & oils

Task

Find out the amount of KOH necessary for the production of FAME by titrating the oil with NaOH stock solution ($c = 1\text{ g/l}$) and separate the oil into its two components glycerol and biodiesel (FAME).


Setup

Titration



- beaker
- test tubes
- 10 ml isopropyl alcohol 99 %
- at least 1 ml of oil (45 °C – for titration)
~10 ml for transesterification
- phenolphthalein solution (~0.3 %)
- water
- NaOH stock solution ($c = 1\text{ g/l}$)
- 2 x one-way syringes with 5 ml
- pipette
- KOH aqueous stock solution (1 g/l)
- bunsen burner/other heat source
- thermometer
- rubber septum / test tube stopper
- balance (only necessary for fresh preparation of the KOH solution)

FAME production





1.1 Production of Biodiesel (FAME) from fats & oils

Performing

Titration:

1. Measure exactly 1 ml of oil with the one-way syringe/pasteur pipette and add approximately 10 ml isopropyl alcohol to the oil into the beaker.
2. Shortly heat the mixture up to 40 °C. Don't bring it to the boil!
3. Add 3 drops of phenolphthalein to the mixture.
4. Now add the NaOH standard solution by a slight swirling until a permanent pink colour appears in the beaker (at least for 30 seconds). Note the amount used.

Used NaOH standard solution in ml: _____

Hint: At this point, all free fatty acids of the sodium hydroxide solution have been converted into its salts and the pH has increased to 8.5 due to the freely existing OH^- ions. The consumption in ml KOH stock solution now has to be converted to the needed amount KOH according to the table/formula. To ensure exact results, the titration may be conducted several times.

As a rule of thumb, the oil quality can be expressed according to its consumption as follows:

- 0 – 3 ml excellent quality; suitable for FAME production
- 3 – 6 ml medium quality; still suitable for FAME production
- > 6 ml low quality; unsuitable for FAME production

5. Using the table, now determine the additional necessary amount of KOH for 1 litre of FAME formation. Determination of the additional necessary amount of NaOH for 1 litre FAME formation.

Consumption of stock solution during titration in ml	Additional necessary amount of KOH in g	Total amount of KOH 5 g + X in g	Acid value mg/g oil (density values for the different oils: cf. exp. 1.3 FAME \leftrightarrow cooking oil, here: 0.92 kg/l)
1.0	1	6	6.5
2.0	2	7	7.6
4.5	4.5	9.5	10.3

6. Fill the necessary amount of methanol and the calculated amount of KOH into the first test tube (TT1) and dissolve it by repeated shaking. It is advisable to calculate the formation for approximately 16 ml of oil due to the typical size of a test tube of 20 ml. Then, enough FAME is produced for all following experiments. The calculation of the proportions is conducted according to the following table:

	Formation for 1000 ml of oil	Formation for 16 ml of oil
oil	1000 ml	16 ml
KOH	5 g + X g (e.g. X = 1 g)	0.1 g + X g
methanol	220 ml	3.5 ml



1.1 Production of Biodiesel (FAME) from fats & oils

Performing

7. Now fill the desired amount of oil for transesterification into test tube 1 (TT1).
8. Heat the mixture in a bain-marie to approx. 50 °C for approx. 20 minutes and shake the test tube from time to time using a stopper.
9. Let the mixture cool down and wait for phase separation (can take up to one hour).

Hint: The heavier brownish glycerol (density 1.2) settles on the ground, while the yellowish biodiesel (density 0.9) swims on top.

10. Wait for phase separation and pipette the swimming biodiesel carefully into a new test tube (TT2). Add water to the biodiesel in TT2 until the test tube is nearly full. Carefully turn the test tube several times by 180° head first to wash away excess methanol and catalyst rests. It is also possible to shake it carefully.

Hint: If the test tube is shaken too heavily, the separation of the two phases may take very long. To detach the solid parts, it can be helpful to additionally filter the FAME/water mixture if the biodiesel does not get clearer during washing.




2.3 Relation between reaction speed and temperature

Task

Investigate the relation between reaction speed and temperature!

Setup



Equipment

- measuring cylinder
- thermometer
- hose
- container for bain-marie

Performing

1. Verify the Q_{10} value with the fermentation experiment by measuring the volume of CO_2 created in a certain time interval once at 20 °C and once at 30 °C (the formation should already be brought to the correct temperature one hour in advance in order to accelerate the reaction). Therefore, fill the provided measuring cylinder with water and put it upside down in a bain-marie so that no air can enter.
2. Then, thread the hose for the proof of CO_2 from below into the measuring cylinder and measure the volume of CO_2 in a certain time interval.
3. From this, calculate the sugar mass that is converted per minute.


Analysis

Temperature in °C	Solubility in g/kg water	mol/kg water	solubility in l/kg water



2.3 Relation between reaction speed and temperature

Analysis



Calculation of sugar translated per minute (in the case of glucose):




4.4 Concentration dependence of the power and the characteristic curve

Task

Check the dependence of the power and the characteristic curve on the concentration.

Setup



Equipment

- motor module
- fuel cell
- 2 hoses
- syringe
- ethanol
- ammeter
- voltmeter

Performing

1. Set up the circuit as given.
2. Dissolve ethanol to 10 % and 20 %.
3. Lead ethanol with different concentrations into the fuel cell. Start with the lower concentrated ethanol, because when starting with higher concentrated ethanol, the fuel cell would saturate and therefore lead to wrong measurement values.
4. Draw both characteristic curves and the voltage-power-diagram into each one diagram and explain the observed trend!

Hint: The cell should be flushed with water before changing the concentration.

Analysis

10 %			20 %		
V in mV	I in mA	P in mW	V in mV	I in mA	P in mW





4.4 Concentration dependence of the power and the characteristic curve

Analysis

10 %			20 %		
V in mV	I in mA	P in mW	V in mV	I in mA	P in mW

Diagrams





4.4 Concentration dependence of the power and the characteristic curve

Analysis

Explanation:

leXsolar GmbH
Strehleener Straße 12-14
01069 Dresden / Germany

Telefon: +49 (0) 351 - 47 96 56 0
Fax: +49 (0) 351 - 47 96 56 - 111
E-Mail: info@lexsolar.de
Web: www.lexsolar.de